

ACE Stir Shaft Quick Reference Chart

Ace Code	Size		Type	Material	Recommended Use																	
	O.D., mm	L, mm			Bearing	Chuck	Agitators															
9534-04	5	318	Button Bottom	Solid Ground Borosilicate Glass	9524-04, 9524-06, 9524-08, 9527-08, 9527-12, 9527-14, 9529	8124, 8126	9542															
9534-06	6	318		Solid Polished Borosilicate Glass																		
9534-40		318		Hollow Borosilicate Glass																		
9535-06	5	318	Knob Bottom	Solid Ground Borosilicate Glass	9524-04, 9527-08	8124, 8126	9541-06															
9541-04		318		PTFE Coated Stainless Steel																		
9541-15		318																				
8134-15	9	380	With Paddles	Solid Polished Borosilicate Glass	8098, 8099, 8133, 9528	8124, 8126	Attached															
8134-25		610	Button Bottom				Solid Ground Borosilicate Glass	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8083, 8085, 8086													
9530	416	Multi-Blade			9530																	
8068-02	10	580	Button Bottom	Solid Ground Borosilicate Glass	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8125, 8126	8082, 8083, 8085, 8086, 8087, 8096															
8068-03		440					Double Multi-Blade	Hollow Ground Borosilicate Glass	Attached													
8068-04		690	Button Bottom	8082, 8083, 8085, 8086, 8087, 8096																		
8068-06		740		Bead Bottom					8088, 8089, 8090													
8068-08		440	Button Bottom				Solid Ground Borosilicate Glass	8082, 8083, 8085, 8086, 8087, 8096														
8068-17		580		Button Bottom				PTFE Coated Glass	8088, 8089, 8090													
8068-18		440	Bottom Drilled Hole				PTFE Coated Stainless Steel		8082, 8083, 8085, 8086, 8087, 8096													
8068-25		440		Button Bottom				Hollow Ground Borosilicate Glass	8088, 8089, 8090													
8068-27		580	Bead Bottom				Solid Ground Borosilicate Glass		8082, 8083, 8085, 8086, 8087, 8096													
8068-30		440		Button Bottom				PTFE Coated Glass	8088, 8089, 8090													
8068-31		690	Bottom Drilled Hole				PTFE Coated Stainless Steel		8082, 8083, 8085, 8086, 8087, 8096													
8068-32		580		Button Bottom				Hollow Ground Borosilicate Glass	8088, 8089, 8090													
8070-05		440	Button Bottom				PTFE Coated Glass		8082, 8083, 8085, 8086, 8087, 8096													
8070-10		690		Bottom Drilled Hole				PTFE Coated Stainless Steel	8088, 8089, 8090													
8071-05		460	Button Bottom				Hollow Ground Borosilicate Glass		8082, 8083, 8085, 8086, 8087, 8096													
8071-07		640		Bottom Drilled Hole				PTFE Coated Stainless Steel	8088, 8089, 8090													
8071-10		690	Button Bottom				Hollow Ground Borosilicate Glass		8082, 8083, 8085, 8086, 8087, 8096													
8073		410		Crescent Shaft				Solid Ground Borosilicate Glass	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 8066, 13443	Attached												
8074-02	420	Remove Button Bottom	Stainless Steel	8044, 8050, 8066, 13443, 13445	8124, 8125, 8126	8082, 8083, 8085, 8086, 8087, 8094, 8095, 8096																
8074-04	450					Button Bottom	Solid Polished Borosilicate Glass	8082, 8083, 8085, 8086, 8087, 8096														
8074-07	580							Plain Shaft	Solid Polished Borosilicate Glass	8094, 8095												
8075-12	440	Button Bottom	Solid Polished Borosilicate Glass	8044, 8066																		
8075-14	580			Plain Shaft	Solid Polished Borosilicate Glass	8044, 8050, 8066, 13443, 13445																
8075-15	690	Bead Bottom	Solid Polished Borosilicate Glass			8044, 8050, 8066, 13443, 13445																
8075-21	440			Button Bottom	Hollow Polished Borosilicate Glass	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 8066, 13443																
8075-23	580	Plain Shaft	Solid Polished Borosilicate Glass			8036, 8038, 8039, 8040, 0, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443																
8075-24	690			Bead Bottom	Solid Polished Borosilicate Glass	Attached																
8075-32	440	Button Bottom	Hollow Polished Borosilicate Glass			8085																
8075-34	580			Bead Bottom	Solid Polished Borosilicate Glass	8091, 8092, 8093, 8100																
8075-36	690	Button Bottom	Hollow Polished Borosilicate Glass			8085																
9523	445			Complete W/Vanes	Hollow Polished Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
9533	475	Button Bottom	Hollow Polished Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-05	700			Button Bottom	Hollow Polished Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-07	900	Bottom Drilled Hole	Solid Polished Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-10	1200			Bottom Drilled Hole	Solid Polished Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-40	900	Bottom Drilled Hole	Solid Polished Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-42	800			Bottom Drilled Hole	Solid Polished Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-43	800	Bottom Drilled Hole	Solid Polished Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-44	1200			Bottom Drilled Hole	Solid Polished Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-45	1140	Bottom Drilled Hole	Solid Polished Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-46	1400			Bottom Drilled Hole	Solid Polished Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8076-48	1400	Bottom Drilled Hole	Solid Polished Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8077-23	700			Button Bottom	Solid Ground Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8077-25	900	Bottom Knob	Hollow Ground Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8077-27	1200			Bottom Knob	Hollow Ground Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8078-05	900	Bottom Knob	Hollow Ground Borosilicate Glass			8050, 8059, 8060, 8061, 8065, 13445	8085															
8078-10	1200			Bottom Knob	Hollow Ground Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8085															
8079-03	700	Bottom Drilled Hole	PTFE Coated Stainless Steel			8059, 8060, 8061, 8065	8085															
8079-05	900			Bottom Drilled Hole	PTFE Coated Stainless Steel	8059, 8060, 8061, 8065	8085															
8079-10	1200	Bottom Drilled Hole	PTFE Coated Stainless Steel			8059, 8060, 8061, 8065	8085															
8080-12	1010			Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101														
8080-14	1140	Bottom Drilled Holes	Solid Polished Borosilicate Glass					8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101												
8080-18	1320									Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101								
8080-22	1470													Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101				
8080-24	1400																	Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101
8080-25	2030																					Bottom Drilled Holes
8080-30	1600			Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157															