

# MOISTURE TEST RECEIVER ...

## *Jacketed*



**New!**

**Improves  
Distillation  
Performance**

### **7745 MOISTURE TEST RECEIVER, Barrett Type, Jacketed ♠**

Used to determine the water content in petroleum and bituminous material. Similar to ASTM Type, except supplied with joints, stopcock to empty measuring tube during the determination and vacuum jacketed side tube to prevent condensation and improve performance. The 10 mL size is subdivided in 0.1 mL from 0-1 mL, 0.2 mL from 1-10 mL; 20 mL size in 0.1 mL from 0-1 mL, 0.2 mL from 1-20 mL. With 2mm bore glass or 1:5 solid PTFE stopcock plug. Joints are \$24/40.



### **Ordering Information:**

Plug Style	Capacity, mL	Ordering Code	Each
Glass	10	<b>7745-102</b>	210.13
PTFE	10	<b>7745-202</b>	226.16
Glass	20	<b>7745-104</b>	216.42
PTFE	20	<b>7745-206</b>	241.81

CONTAINS  
2010  
PRICING

♠ACE manufactured item

INTERNATIONAL  
SALES

Phone: 856.692.3333  
Fax: 856.692.8919

e-mail:  
export@aceglass.com



**ACE GLASS INCORPORATED**

P.O. Box 688 • Vineland, NJ 08362-0688 • 856-692-3333 • Fax: 856-692-8919

**TOLL-FREE: 1-800-223-4524 • FAX: 1-800-543-6752**

**www.aceglass.com email: sales@aceglass.com**

11/2/2009

Prices and specifications are  
subject to change without notice