

# BIG JARS WITH HANDLES, AND GRADUATED

**NEW!**

RECESSED HANDLE



RECESSED HANDLE




## 6233 JARS Cylindrical with Handles ★

Mold-blown, heavy walled jars with open ground flat top and recessed finger grip handles for easier handling. Graduated. Not recommended for use with direct heat.

Capacity, Liters	O.D., mm	Height, mm	Order Code	Each
7.25	223 (8.75")	254 (10")	<b>6233-07</b>	436.04
9.25	257 (10.13")	254 (10")	<b>6233-09</b>	498.35
17.0	305 (12")	305 (12")	<b>6233-17</b>	568.57
26.5	305 (12")	457 (18")	<b>6233-26</b>	696.55
32.0	406 (16")	305 (12")	<b>6233-32</b>	1,130.34

6228,6231



## 6228, 6231 BEAKER Pilot Plant

Large sizes for batch operations and mixing large volumes of measured liquids.

Capacity, mL	Subdiv., mL	Approx. Height, mm	O.D., mm	Order Code	Each
5,000	100	457	152	<b>6228-05</b>	▲297.96
10,000	100	457	223	<b>6228-10</b>	▲471.64
15,000	500	390	260	<b>6231-21</b>	★444.79
20,000	500	430	285	<b>6231-27</b>	★489.39



**ACE GLASS INCORPORATED**

P.O. Box 688 • Vineland, NJ 08362-0688 • 856-692-3333 • Fax: 856-692-8919

TOLL-FREE: 1-800-223-4524 • FAX: 1-800-543-6752

www.aceglass.com email: sales@aceglass.com

11/2/2009

Prices and specifications are subject to change without notice

★ Net