





[www.DelpMixer.com](http://www.DelpMixer.com)  
[Info@DelpMixer.com](mailto:Info@DelpMixer.com)  
1-800-504-1935

2528 Billingsley Road  
Columbus, Ohio 43235  
614-889-7331 (Fax)

## 316 Series Mixers



**316S**



**316AF**

- Equipped with a 0.50 in. diameter shaft and an extra bearing for heavy duty mixing
- Operates in pressure (200psi max) or vacuum (30 in. Hg) environments
- Manufactured from high-grade materials to withstand autoclave sterilization
- 316S with variable speed electric motor provides 22 in/oz (1.375 in/lb or 1.58 kg/cm) of torque
- 316AF for standard cable or direct-drive electric, air or hydraulic drive motors; capable of up to 132 in/oz (8.25 in/lb or 9.50 kg/cm) of torque with standard rare earth magnets