

JACKETED *EUROPEAN STYLE/TAPERED* FLASKS

More Efficient Stirring while Heating/Cooling

6959 FLASK *European Style, Three Necks, Jacketed* ★

Fully jacketed European style, tapered walls with shallow bottom for more efficient magnetic or mechanical half-moon shaped paddle stirring. Easy to clean. Supplied with ₤ center neck, two ₤ side necks.

Capacity, mL	₤ Center Neck	₤ Side Necks	Size I.D. Tube for H/C, mm	Order Code	Each
125	14/20	14/20	9.5" (3/8")	6959-03	297.36
125	24/40	14/20	9.5" (3/8")	6959-05	297.36
250	14/20	14/20	11.1" (7 ³ / ₁₆ ")	6959-11	333.77
250	24/40	24/40	11.1" (7 ³ / ₁₆ ")	6959-13	333.77
500	14/20	14/20	12.7 (1/2")	6959-22	400.53
500	24/40	24/40	12.7 (1/2")	6959-24	400.53
500	29/42	24/40	12.7 (1/2")	6959-26	400.53
1000	24/40	24/40	12.7 (1/2")	6959-30	423.59
1000	29/42	24/40	12.7 (1/2")	6959-33	423.59
2000	24/40	24/40	12.7 (1/2")	6959-36	603.21
2000	29/42	24/40	12.7 (1/2")	6959-38	603.21



6959-13

6959 FLASK *European Style, Five Necks, Jacketed* ★

Fully jacketed European style, tapered walls with shallow bottom for more efficient magnetic or mechanical half-moon shaped paddle stirring. Easy to clean. Supplied with ₤ center neck, two ₤ side necks, and two Ace-Threds for thermometers, etc., up to 7 mm O.D. Supplied with nylon bushings and FETFE O-Rings.

Capacity, mL	₤ Center Neck	₤ Side Necks	Size I.D. Tube for H/C, mm	Order Code	Each
125	14/20	₤14/20-#7	9.5" (3/8")	6959-44	361.67
250	24/40	₤14/20-#7	11.1" (7 ³ / ₁₆ ")	6959-48	461.20
250	29/42	₤14/20-#7	11.1" (7 ³ / ₁₆ ")	6959-51	461.20
500	24/40	₤24/40-#7	12.7 (1/2")	6959-55	515.83



6959-48

6960 FLASK *Tapered Wall, Five Necks, Jacketed* ★

Fully jacketed with flask walls tapered toward the bottom. This design allows for continuous operation with smaller volumes. Supplied with ₤ center neck, two ₤ side necks, and two #7 Ace-Threds for thermometers, etc., up to 7 mm O.D. Supplied with nylon bushings and FETFE O-Rings.

Capacity, mL	₤ Center Neck	₤ Side Necks	Size I.D. Tube for H/C, mm	Order Code	Each
50	14/20	₤14/20-#7	9.5" (3/8")	6960-02	423.59
125	14/20	₤14/20-#7	9.5" (3/8")	6960-09	423.59
250	24/40	₤24/40-#7	11.1" (7 ³ / ₁₆ ")	6960-18	479.41
250	29/42	₤24/40-#7	11.1" (7 ³ / ₁₆ ")	6960-20	479.41
500	24/40	₤24/40-#7	12.7 (1/2")	6960-35	576.51
500	29/42	₤24/40-#7	12.7 (1/2")	6960-37	576.51



6960-35



P.O. Box 688 • Vineland, NJ 08362-0688 • 856-692-3333 • Fax: 856-692-8919

★Net

TOLL-FREE: 1-800-223-4524 • FAX: 1-800-543-6752

www.aceglass.com email: sales@aceglass.com

01/6/2011

Prices and specifications are subject to change without notice