




ILMVAC

**New!**

INTRODUCES...

# Diaphragm Vacuum Pumps

*Ace Glass is pleased to announce the addition of four new oil-free, chemically resistant diaphragm vacuum pumps. Ideal for Rotary Evaporators and other acid or solvent vapor applications. The pumps are compact and extremely quiet — great for all laboratory applications. Low maintenance due to beltless, oil-free operation, and all wetted parts are PTFE or PTFE-type compounds. This selection of pumps features wide flow rate and vacuum ranges, and all come with a two-year manufacturer's warranty.*



14112-07

## MPC101Z

### 14112-07 Vacuum Pump ★

Two-stage, chemically resistant diaphragm type pump. Ultimate pressure 6 torr (8 mBar), 18L/min. flow rate. Connections are for 8mm I.D. vacuum hose. 115v, 50/60hz.

Order Code	Each
14112-07	1,739.44

## MPC104T

### 14112-09 Vacuum Pump ★

Three-stage, chemically resistant diaphragm type pump. Ultimate pressure 1.5 torr (2 mBar), 16L/min. flow rate. Connections are for 8mm I.D. vacuum hose. 115v, 50/60hz.

Order Code	Each
14112-09	2,164.90



14112-09