



INTRODUCES...

# Suba-Septa™

*Suba-Seal™ Septa, the highest quality roll-over style septa in the world. The bottom of the septa is serrated for maximum sealing even in open-top vessels and containers. Certain sizes are also perfect for maximum sealing in standard taper joints. These septa are ideal for moisture or air sensitive applications. The septa when pierced with a non-coring needle have excellent re-sealing capability. The septa have double sealing feature as the turn-over top seals on the outside of the vessel and the serrated bottom flanges seal excellently to the inner surface of the container. Suba-Seal septa are autoclaveable and are packaged in “white room” conditions.*



## 9107 SUBA-SEAL™ SEPTA Red Rubber ★

Description	Fits	Pack Size	Order Code	Each
Sample kit, 10 of each size	All	Sample Set	<b>9107-04</b>	208.35
#9 fits neck ID 8mm		100	<b>9107-116</b>	59.99
#13 fits neck ID 9.5mm	10/30	100	<b>9107-110</b>	61.39
#17 fits neck ID 11mm		100	<b>9107-114</b>	69.36
#21 fits neck ID 12.5mm		100	<b>9107-118</b>	74.82
#29 fits neck ID 16mm		100	<b>9107-126</b>	112.15
#33 fits neck ID 17.5mm		100	<b>9107-130</b>	121.51
#37 fits neck ID 19mm	19/22	100	<b>9107-134</b>	127.26
#41 fits neck ID 20.5mm		100	<b>9107-138</b>	155.22
#45 fits neck ID 22mm	24/40	100	<b>9107-142</b>	176.19
#49 fits neck ID 24mm		100	<b>9107-146</b>	224.44
#53 fits neck ID 25.5mm		100	<b>9107-150</b>	251.70
#57 fits neck ID 27mm	29/42	100	<b>9107-154</b>	290.86



★ Prices are NET.

Suba-Seal™ is a registered trademark of Sigma-Aldrich Biotechnology LP.