

O-RING KITS

Viton® A and Buna-N

Ace Glass is pleased to offer O-Ring kits in Buna-N, and Viton. The O-Rings are packaged in a clear, plastic box with a resealable lid, clearly marked with material type, O-Ring size, and quantity. The cardboard insert indicates the individual Ace codes for each size, so you can easily re-order the O-Rings as they are depleted. The quantity of each O-Ring varies by size and compartment.



O-Ring Size	Dimensions (mm)	
	I.D.	Width
-006	2.9	1.78
-007	3.7	1.78
-008	4.5	1.78
-009	5.3	1.78
-010	6.1	1.78
-011	7.7	1.78
-012	9.3	1.78
-110	9.2	2.6
-111	10.8	2.6
-112	12.4	2.6
-113	13.9	2.6
-114	15.5	2.6
-115	17.1	2.6
-116	18.7	2.6
-210	18.6	3.5
-211	20.2	3.5
-212	21.8	3.5
-213	23.4	3.5
-214	25.0	3.5
-215	26.6	3.5
-216	28.2	3.5
-217	29.7	3.5
-218	31.3	3.5
-219	32.9	3.5
-220	34.5	3.5
-221	36.1	3.5
-222	37.7	3.5
-325	37.5	5.3
-326	40.6	5.3
-327	43.8	5.3

Viton® A Kit 7855-99

- Comprised of linear co-polymer of vinylidene fluoride and hexafluoropropylene
- Fluorine content is between 66 to 70%
- Generally black in color, and with generally good temperature and chemical resistance
- Kit contains 500 total o-rings

Buna-N Kit 7855-499

- Also known as Nitrile rubber
- A synthetic blended rubber, co-polymer of acrylonitrile and butadiene
- Generally resistant to fuels and oils
- Generally black or brown in color
- Kit contains 500 total o-rings

Viton® A and Buna-N Kit contents

-006 Qty. 35	-007 Qty. 35	-008 Qty. 35	-009 Qty. 35	-010 Qty. 35	-011 Qty. 35	-012 Qty. 35	-110 Qty. 20	-111 Qty. 20
-112 Qty. 20	-113 Qty. 20	-114 Qty. 16	-115 Qty. 16	-116 Qty. 16	-210 Qty. 10	-211 Qty. 10	-212 Qty. 10	-213 Qty. 10
-214 Qty. 10	-215 Qty. 10	-216 Qty. 10	-217 Qty. 7	-218 Qty. 7	-219 Qty. 7	-220 Qty. 7	-221 Qty. 7	-222 Qty. 7
-325 Qty. 5		-327 Qty. 5				-326 Qty. 5		

FETFE & Mixed Material O-Ring Kits

Ace Glass is pleased to offer O-Ring kits in FETFE, and a mixed set (FETFE, Silicone, Buna-N and Viton A). The O-Rings are packaged in a clear, plastic box with a resealable lid, clearly marked with material type, O-Ring size, and quantity. The cardboard insert indicates the individual Ace codes for each size, so you can easily re-order the O-Rings as they are depleted. The quantity of each O-Ring varies by size and compartment.



O-Ring Size	Dimensions (mm)	
	I.D.	Width
-008	4.5	1.78
-009	5.3	1.78
-010	6.1	1.78
-011	7.7	1.78
-012	9.3	1.78
-013	10.8	1.78
-014	12.4	1.78
-018	18.8	1.78
-108	6.0	2.6
-110	9.2	2.6
-111	10.8	2.6
-112	12.4	2.6
-116	18.7	2.6
-121	26.6	2.6
-136	50.5	2.6
-210	18.6	3.5
-212	21.8	3.5
-220	34.5	3.5

FETFE® Kit 7855-799

- An ACE exclusive compound which is a fluoroelastomer with TFE additives
- Good compression set along with excellent temperature and chemical compatibility
- Generally black or white in color
- Kit contains 216 total o-rings, 12 of each size

-008	-010	-011	-009	-012	-013
-108	-014	-018	-110	-111	-112
-116	-121	-136	-210	-212	-220

Mixed Kit 7855-300

- Includes 4 different materials: FETFE, Silicone, Buna-N, Viton A
- Kit contains 360 total o-rings

Viton A -010	Viton A -011	Viton A -013	Viton A -108	Viton A -110	Silicone -010
Silicone -011	Silicone -013	Silicone -110	Silicone -210	FETFE -008	FETFE -013
FETFE -110	FETFE -210	FETFE -212	Buna-N -011	Buna-N -013	Buna-N -110

