

ACE Stir Shaft Quick Reference Chart

Ace Code	Size		Type	Material	Recommended Use																																	
	O.D., mm	L, mm			Bearing	Chuck	Agitators																															
9534-04	5	318	Button Bottom	Solid Ground Borosilicate Glass	9524-04, 9524-06, 9524-08, 9527-08, 9527-12, 9527-14, 9529	8124, 8126	9542																															
9534-06	6	318		Solid Polished Borosilicate Glass																																		
9534-40	6	318		Hollow Borosilicate Glass																																		
9535-06		318	Knob Bottom	Solid Ground Borosilicate Glass	9524-04, 9527-08		8124, 8126	9541-06																														
9541-04	5	318		PTFE Coated Stainless Steel				9541-07																														
9541-15	5	318						Attached																														
8134-15		380	With Paddles					8083, 8085, 8086																														
8134-25	9	610	Button Bottom	Solid Polished Borosilicate Glass	8098, 8099, 8133, 9528		8124, 8126	9530																														
9530-04		416		Solid Ground Borosilicate Glass				8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	8082, 8083, 8085, 8086, 8087, 8096																												
8068-02		580								Multi-Blade	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	Attached																									
8068-03		440				Double Multi-Blade							8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	8082, 8083, 8085, 8086, 8087, 8096																							
8068-04		690								Button Bottom	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126			8088, 8089, 8090, 8091																							
8068-06		740				Bead Bottom							8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	8082, 8083, 8085, 8086, 8087, 8096																							
8068-08		440								Button Bottom	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126			8088, 8089, 8090, 8091																							
8068-17		580		Bead Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8082, 8083, 8085, 8086, 8087, 8096																													
8068-18		440							Button Bottom	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	8088, 8089, 8090, 8091																										
8068-25		440		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126				8082, 8083, 8085, 8086, 8087, 8096																										
8068-27		580	Crescent Shaft		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8088, 8089, 8090, 8091																													
8068-30		440		Remove Button Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	Attached																													
8068-31		690	Button Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8082, 8083, 8085, 8086, 8087, 8096																													
8068-32		580		Plain Shaft		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8094, 8095																													
8070-05		440	Bead Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8088, 8089, 8090, 8091																													
8070-10		690		Complete W/Vanes		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	Attached																													
8071-05		460	Hollow Polished Borosilicate Glass		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8091, 8092, 8093, 8100																													
8071-07		640		Solid Polished Borosilicate Glass		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8071-10		690	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8091, 8092, 8093, 8100																													
8073	10	410		Button Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8074-02		420	Plain Shaft		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8091, 8092, 8093, 8100																													
8074-04		460		Bead Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8074-07		580	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8075-12		440		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8075-14		580	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8075-15		690		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8075-21		440	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8075-23		580		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8075-24		690	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8075-32		440		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8075-33		500	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8075-34		580		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8075-36		690	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8075-38		560		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
9523-04		445	Complete W/Vanes		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
9533-02		475		Hollow Polished Borosilicate Glass		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8076-03	19	560	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8076-05		710		Hollow Polished Borosilicate Glass		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8076-07		910	Button Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8076-10		1210		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8076-39		940	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8076-40		910		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8076-41		990	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8076-42		810		Bottom Knob		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8076-43		810	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8076-44		1210		Bottom Knob		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8076-45		1140	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8076-46		1380		Bottom Knob		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8076-48		1380	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8077-23		710		Button Bottom		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8077-25		910	Bottom Knob		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8077-27		1210		Bottom Knob		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8078-05		900	Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																													
8078-10		1200		Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126	8085																													
8079-03	730	Bottom Drilled Hole	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126		8085																															
8079-05	910			Bottom Drilled Hole		8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	8085																														
8079-10	1210	Bottom Drilled Holes	8036, 8038, 8039, 8040, 8041, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443		8124, 8126			8085																														
8080-12	1010			Bottom Drilled Holes		Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101																													
8080-14	1140	Bottom Drilled Holes	Solid Polished Borosilicate Glass		8050, 8067				6472-155, 6472-156, 6472-157	8091, 8093, 8101																												
8080-16	1300									Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101																								
8080-18	1320													Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101																				
8080-22	1470																	Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101																
8080-24	1400																					Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101												
8080-25	2080																									Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101								
8080-29	1900																													Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101				
8080-30	1600																																	Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	6472-155, 6472-156, 6472-157	8091, 8093, 8101