

## ACE Stir Shaft Quick Reference Chart

Ace Code	Size		Type	Material	Recommended Use			
	O.D., mm	L, mm			Bearing	Chuck	Agitators	
9534-04	5	318	Button Bottom	Solid Ground Borosilicate Glass	9524-04, 9524-06, 9524-08, 9527-08, 9527-12, 9527-14, 9529	8124, 8126	9542	
9534-06	6	318		Solid Polished Borosilicate Glass				
9534-40	6	318		Hollow Borosilicate Glass				
9535-06		318	Knob Bottom	Solid Ground Borosilicate Glass	9524-04, 9527-08		8124, 8126	9541-06
9541-04	5	318		PTFE Coated Stainless Steel				9541-07
9541-15	5	318						Attached
8134-15		380	With Paddles					8083, 8085, 8086
8134-25	9	610	Button Bottom	Solid Polished Borosilicate Glass	8098, 8099, 8133, 9528		8124, 8126	9530
9530-04		416						
8068-02		580						
8068-03		440						
8068-04		690						
8068-06		740						
8068-08		440		Multi-Blade				
8068-17		580		Double Multi-Blade				
8068-18		440						
8068-25		440		Button Bottom		Hollow Ground Borosilicate Glass		8036, 8038, 8039, 8040, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443
8068-27		580						
8068-30		440	Bead Bottom	Solid Ground Borosilicate Glass	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8051, 8053, 8055, 13443	8124, 8126	Attached	
8068-31		690						
8068-32		580	Button Bottom	PTFE Coated Glass	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8051, 8053, 8055, 13443	8124, 8126	8082, 8083, 8085, 8086, 8087, 8096	
8070-05		440						
8070-10		690	Bottom Drilled Hole	PTFE Coated Stainless Steel	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8051, 8053, 8055, 13443	8124, 8126	8088, 8089, 8090, 8091	
8071-05		460						
8071-07		640						
8071-10		690						
8073	10	410	Crescent Shaft	Solid Ground Borosilicate Glass	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	Attached	
8074-02		420	Remove Button Bottom	Stainless Steel	8044, 8050, 8066, 13443, 13445	8124, 8126	8082, 8083, 8085, 8086, 8087, 8094, 8095, 8096	
8074-04		460						
8074-07		580	Button Bottom	Solid Polished Borosilicate Glass	8044, 8066	8124, 8126	8082, 8083, 8085, 8086, 8087, 8096	
8075-12		440						
8075-14		580	Plain Shaft	Solid Polished Borosilicate Glass	8044, 8066	8124, 8126	8094, 8095	
8075-15		690						
8075-21		440	Bead Bottom	Solid Polished Borosilicate Glass	8044, 8050, 8066, 13443, 13445	8124, 8126	8088, 8089, 8090, 8091	
8075-23		580						
8075-24		690	Complete W/Vanes	Hollow Polished Borosilicate Glass	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 8066, 13443	8124, 8126	Attached	
8075-32		440						
8075-33		500	Complete W/Vanes	Solid Ground Borosilicate Glass	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 13443	8124, 8126	Attached	
8075-34		580						
8075-36		690	Bottom Drilled Hole	Solid Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8091, 8092, 8093, 8100	
8075-38		560						
9523-04		445	Complete W/Vanes	Hollow Polished Borosilicate Glass	8036, 8038, 8039, 8040, 8042, 8043, 8047, 8050, 8051, 8053, 8055, 8066, 13443	8124, 8126	Attached	
9533-02		475						
8076-03	19	560	Bottom Drilled Hole	Solid Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8091, 8092, 8093, 8100	
8076-05		710						
8076-07		910	Button Bottom	Hollow Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8085	
8076-10		1210						
8076-39		940	Bottom Drilled Hole	Solid Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8091, 8092, 8093, 8100	
8076-40		910						
8076-41		990	Bottom Drilled Hole	Solid Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8091, 8092, 8093, 8100	
8076-42		810						
8076-43		810	Bottom Drilled Hole	Solid Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8091, 8092, 8093, 8100	
8076-44		1210						
8076-45		1140	Bottom Drilled Hole	Solid Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8091, 8092, 8093, 8100	
8076-46		1380						
8076-48		1380	Bottom Drilled Hole	Solid Polished Borosilicate Glass	8049, 8050, 8067, 13445	8124, 8126	8091, 8092, 8093, 8100	
8077-23		710						
8077-25		910	Button Bottom	Solid Ground Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8124, 8126	8085	
8077-27		1210						
8078-05		900	Bottom Knob	Hollow Ground Borosilicate Glass	8050, 8059, 8060, 8061, 8065, 13445	8124, 8126	8091, 8092, 8093, 8100	
8078-10		1200						
8079-03		730	Bottom Drilled Hole	PTFE Coated Stainless Steel	8059, 8060, 8061, 8065	8124, 8126	8091, 8092, 8093, 8100	
8079-05		910						
8079-10		1210	Bottom Drilled Holes	Solid Polished Borosilicate Glass	8050, 8067	8124, 8126	8091, 8093, 8101	
8080-12	28	1010						
8080-14		1140						
8080-16		1300						
8080-18		1320						
8080-22		1470						
8080-24		1400						
8080-25		2080						
8080-29		1900						
8080-30		1600						