
IMPORTANT INFORMATION

ABOUT 5555 STORAGE BOTTLE

Please Read Before Using

5555 STORAGE BOTTLE, Heavy Wall, Safety Coated

These bottles are intended only for storage of low-boiling liquids at room temperature or at not less than -20°C . They are designed to resist internal pressures up to 200 psig, as low as -20°C . The low temperature rating is based primarily on the rating for PTFE.

Due to varying conditions, ACE cannot guarantee these vessels from breakage under pressure. **ALL LABORATORY SAFETY PROCEDURES SHOULD BE OBSERVED.** Always work behind a shield.

Safety coating helps prevent shattering and reduces spills. Although coating reduces chances of scratching, bottles should not be stacked on top of each other nor touching each other when subjected to vibration or shock.

Do not remove frost by placing in hot water, etc.

Do not use with materials which solidify on standing and create excessive stress on the glass.

***ALWAYS OBSERVE PROPER
LABORATORY SAFETY PROCEDURES***



ACE GLASS INCORPORATED

P.O. BOX 688, VINELAND, NJ 08380-0688 • 800-223-1624 • FAX 800-543-8762