

**Table 1**

<b>Meter Tube No.</b>	<b>Float</b>	<b>Flow Range mL/min.</b>	
		<b>Air Flow</b>	<b>Water Flow</b>
31	Stainless Steel Glass	5.0-45 None	0.065-0.65 —
32	Stainless Steel Glass	120-1200 48-475	03.0-30 0.5-6.0
33	Stainless Steel Glass	475-5700 240-2800	10-180 5-60
34	Stainless Steel Glass	950-17,000 None	50-500 —
35	Stainless Steel Glass	4.8-50L/min. 2.5-28L/min.	— —
36 Wide Range	Stainless Steel  Glass	None 0.24-14.3 L/min.	— 6.5-300