



# 100mL to 6L BENCH-TOP STANDS

Stand only, designed for Ace's Scale-Up Series Reactor product line. Stand features: 100mL to 6000mL capacity, lightweight aluminum, adjustable leveling feet and quick adjustment rod-holding brackets for up to 5/8" (16mm) diameter rods. *Not included: our 6442 family mounting clamp and our 11177 family bottom ring support which must be sized to the supported vessel.*

Fits both filter and pressure reactor styles as well as standard round bottom flasks. It will accommodate flasks with or without bottom valves. The motor mount raises and lowers as does the flask support system. This aids in aligning the flask, motor and stir rod. Motor mount accepts overhead stir motors.

## Scale-Up Series Stand, Single

- Compatible with jacketed and unjacketed flasks
- Adjustable leveling feet
- (3) Quick adjusting rod-holding brackets for up to 5/8" O.D. rods

### 100mL to 2L Capacity

Material (Frame)	Aluminum
Dimensions (in)	38 x 24.75 x 19.5 (H x W x D)
Order Code	12841-02


### 3L to 6L Capacity

Material (Frame)	Aluminum
Dimensions (in)	48 x 24.75 x 19.5 (H x W x D)
Order Code	12841-01

### Accessories

Polyethylene Spill Tray - 4.5L Capacity	
Order Code	12841-50

Vertical Mounting Rod Assembly, 36" height	
Order Code	12841-12



## Scale-Up Series Stand, Dual

- Compatible with jacketed and unjacketed flasks
- Adjustable leveling feet
- (6) Quick adjusting rod-holding brackets for up to 5/8" O.D. rods

### 100mL to 2L Capacity

Material (Frame)	Aluminum
Dimensions (in)	38 x 30 x 19.5 (H x W x D)
Order Code	12843-38

### 3L to 6L Capacity

Material (Frame)	Aluminum
Dimensions (in)	48 x 30 x 19.5 (H x W x D)
Order Code	12843-48

### Accessories

Polyethylene Spill Tray - 6L Capacity	
Order Code	12843-50



100mL to 6L

# BENCH-TOP STANDS

Caframo® H-Base stands offer improved stability for benchtop reaction assemblies. The base is manufactured from cast zinc aluminum and coated with a chemical resistant epoxy paint. Each upright support tube is manufactured from 304 stainless steel and plastic capped.

## Caframo® H-Base Stand, Improved Stability

- Compatible with jacketed and unjacketed flasks
- Adjustable leveling feet
- (2) threaded holes for additional 5/8" O.D. support tubes
- Accepts rod mounted motors with 13568-16 Clamp (not included)

### 28in Height

Material (Base)	Cast zinc-aluminum
Material (Upright)	304 stainless steel
Dimensions (in)	28 x 17 x 17 (H x W x D)
Order Code	13568-02

### 38in Height

Material (Base)	Cast zinc-aluminum
Material (Upright)	304 stainless steel
Dimensions (in)	38 x 17 x 17 (H x W x D)
Order Code	13568-04

### 48in Height

Material (Base)	Cast zinc-aluminum
Material (Upright)	304 stainless steel
Dimensions (in)	48 x 17 x 17 (H x W x D)
Order Code	13568-06

### Accessories

Clamp, supports up to 5/8" rod	13568-16
Secondary Support, 28" height	13568-25
Support Rod, 36" x 5/8" O.D.	13586-27



## U-Shaped Stand, Heavy Duty

- Compatible with jacketed and unjacketed flasks
- Can accommodate vessel up to 18" O.D.
- (2) threaded holes for additional 5/8" O.D. support tubes

### 28in Height

Material	303 Stainless steel
Dimensions (in)	28 x 17.75 x 10.5 (H x W x D)
Order Code	13586-10

### 36in Height

Material	303 Stainless steel
Dimensions (in)	36 x 17.75 x 10.5 (H x W x D)
Order Code	13586-13

### Accessories

Support Rod, 60" x 5/8" O.D.	13586-15
Support Rod, 28" x 5/8" O.D.	13586-25
Support Rod, 36" x 5/8" O.D.	13586-27



1430 Northwest Boulevard  
Vineland, NJ 08360

sales@aceglass.com  
P: (800) 223-4524  
F: (800) 543-6752