

COMPRESSION ADAPTERS

Spherical Joint Adapters

These adapters feature a ground glass tube sidearm which allows for the attachment of compression fittings. They offer a versatile connection for adapting flexible or rigid tubing. They may also be used for support of probes or other rigid body equipment. The sidearm is used to make tube compression fitting connection to a wide variety of flexible and rigid tubing using our 12709, 12711, 12715 & 12716 product families of fittings.



**Tube Compression,
Right Angle**

Joint Type	Tube O.D. (in.)	Order Code
28/15 Joint		
Socket	3/8	12719-28
Ball	3/8	12719-30
Socket	1/2	12719-32
Ball	1/2	12719-34
Socket	3/4	12719-36
Ball	3/4	12719-38
35/25 Joint		
Socket	3/8	12719-40
Ball	3/8	12719-42
Socket	1/2	12719-44
Ball	1/2	12719-46
Socket	3/4	12719-48
Ball	3/4	12719-50
DN25 Joint		
Socket	3/8	12719-52
Ball	3/8	12719-54
Socket	1/2	12719-56
Ball	1/2	12719-58
Socket	3/4	12719-60
Ball	3/4	12719-62
DN40 Joint		
Socket	3/8	12719-64
Ball	3/8	12719-66
Socket	1/2	12719-68
Ball	1/2	12719-70
Socket	3/4	12719-72
Ball	3/4	12719-74



**Tube Compression,
Vertical**

Joint Type	Tube O.D. (in.)	Order Code
28/15 Joint		
Socket	3/8	12731-28
Ball	3/8	12731-30
Socket	1/2	12731-32
Ball	1/2	12731-34
Socket	3/4	12731-36
Ball	3/4	12731-38
35/25 Joint		
Socket	3/8	12731-40
Ball	3/8	12731-42
Socket	1/2	12731-44
Ball	1/2	12731-46
Socket	3/4	12731-48
Ball	3/4	12731-50
DN25 Joint		
Socket	3/8	12731-52
Ball	3/8	12731-54
Socket	1/2	12731-56
Ball	1/2	12731-58
Socket	3/4	12731-60
Ball	3/4	12731-62
DN40 Joint		
Socket	3/8	12731-64
Ball	3/8	12731-66
Socket	1/2	12731-68
Ball	1/2	12731-70
Socket	3/4	12731-72
Ball	3/4	12731-74



**Tube Compression,
Twin Right Angle**

Joint Type	Tube O.D. (in.)	Order Code
28/15 Joint		
Socket	3/8	12737-28
Ball	3/8	12737-30
Socket	1/2	12737-32
Ball	1/2	12737-34
Socket	3/4	12737-36
Ball	3/4	12737-38
35/25 Joint		
Socket	3/8	12737-40
Ball	3/8	12737-42
Socket	1/2	12737-44
Ball	1/2	12737-46
Socket	3/4	12737-48
Ball	3/4	12737-50
DN25 Joint		
Socket	3/8	12737-52
Ball	3/8	12737-54
Socket	1/2	12737-56
Ball	1/2	12737-58
Socket	3/4	12737-60
Ball	3/4	12737-62
DN40 Joint		
Socket	3/8	12737-64
Ball	3/8	12737-66
Socket	1/2	12737-68
Ball	1/2	12737-70
Socket	3/4	12737-72
Ball	3/4	12737-74

Standard Taper Joint Compression Adapters

These adapters feature a ground glass tube sidearm which allows for the attachment of compression fittings. They offer a versatile connection for adapting flexible or rigid tubing. They may also be used for support of probes or other rigid body equipment. The sidearm is used to make tube compression fitting connection to a wide variety of flexible and rigid tubing using our 12709, 12711, 12715 & 12716 product families of fittings.




**Tube Compression,
Right Angle**

Joint Size	Tube O.D. (in.)	Order Code	
w/o Drip Tip			
14/20	1/4	12719-02	
	3/8	12719-04	
24/40	1/4	12719-06	
	3/8	12719-08	
	1/2	12719-10	
29/42	3/4	12719-12	
	1/4	12719-14	
	3/8	12719-16	
45/50	1/2	12719-18	
	3/4	12719-20	
	3/8	12719-22	
	1/2	12719-24	
45/50	3/4	12719-26	
	w/Drip Tip		
	14/20	1/4	12722-01
3/8		12722-03	
24/40	1/4	12722-07	
	3/8	12722-09	
	1/2	12722-11	
29/42	3/4	12722-13	
	1/4	12722-15	
	3/8	12722-17	
45/50	1/2	12722-19	
	3/4	12722-21	
	3/8	12722-23	
	1/2	12722-25	
45/50	3/4	12722-27	



**Tube Compression,
Vertical**

Joint Size	Tube O.D. (in.)	Order Code
w/o Drip Tip		
14/20	1/4	12731-02
	3/8	12731-04
24/40	1/4	12731-06
	3/8	12731-08
	1/2	12731-10
29/42	3/4	12731-12
	1/4	12731-14
	3/8	12731-16
45/50	1/2	12731-18
	3/4	12731-20
	3/8	12731-22
45/50	1/2	12731-24
	3/4	12731-26
	w/Drip Tip	
14/20	1/4	12736-01
	3/8	12736-03
24/40	1/4	12736-07
	3/8	12736-09
	1/2	12736-11
29/42	3/4	12736-13
	1/4	12736-15
	3/8	12736-17
45/50	1/2	12736-19
	3/4	12736-21
	3/8	12736-23
45/50	1/2	12736-25
	3/4	12736-27



**Tube Compression,
Twin Right Angle**

Joint Size	Tube O.D. (in.)	Order Code
w/o Drip Tip		
14/20	1/4	12737-02
	3/8	12737-04
24/40	1/4	12737-06
	3/8	12737-08
	1/2	12737-10
29/42	3/4	12737-12
	1/4	12737-14
	3/8	12737-16
45/50	1/2	12737-18
	3/4	12737-20
	3/8	12737-22
45/50	1/2	12737-24
	3/4	12737-26
	w/Drip Tip	
14/20	1/4	12739-01
	3/8	12739-03
24/40	1/4	12739-05
	3/8	12739-07
	1/2	12739-09
29/42	3/4	12739-11
	1/4	12739-13
	3/8	12739-15
45/50	1/2	12739-17
	3/4	12739-19
	3/8	12739-21
45/50	1/2	12739-23
	3/4	12739-25