



SERIES 12842

SUPPORT STAND


Universal open frame stand design allows for maximum clearance and access to the ports on the head. To allow for vessel diversity and future scale-up, the frame has been designed to be fitted with Ace Glass large scale reactors.

Kilo Scale Support Stands:

- 10L to 150L
- 82.25" or 96.25" heights available
- Flange or Rod motor mount
- Flange clamp for vessel
- Adjustable mounting bars
- Support shelf for flask
- Easy access for materials collection
- Optional accessories for large scale distillation

Reactor Size, L	Flange Size, mm	Motor Mount Type	D x W x H in	Order Code
10-20	200	Rod	27.5 x 24.25 x 82.25	12842-02
10-20	200	Flange	27.5 x 24.25 x 82.25	12842-04
30-50	300	Flange	27.5 x 24.25 x 82.25	12842-06
50-150	300	Flange	27.5 x 24.25 x 96.25	12842-08
100 (<i>low profile</i>)	400	Flange	27.5 x 24.25 x 82.25	12842-10
50-150	400	Flange	27.5 x 24.25 x 96.25	12842-12
10-50	N/A	Dual Rod	27.5 x 24.25 x 82.25	12842-20

KILO-SCALE



Also Available...

Contact us for custom large and bench scale applications.




Contact your local sales rep today for more information

1-800-223-4524 | www.aceglass.com

SERIES 12842

SUPPORT ACCESSORIES



Accessory

	Order Code
Bracing Set for Flange Mount, strength for weight of hazardous duty motors	12842-41
Safety Shield (Lexan™ and aluminum), fits 10L to 150L reactors	12842-42
Shelf with (2) 28" x 5/8" Stainless Steel Supports	12842-43
Rear Panel Safety Shield, (Lexan™ and aluminum)	12842-44
Swinging Motor Mount, designed for use with our 13553 standard motor and the (13569-01) Caframo 1540 Crossover	12842-45
Bolt Latch Clamp, 200mm Flange	12842-50
Bolt Latch Clamp, 300mm Flange	12842-60
Bolt Latch Clamp, 400mm Flange	12842-70
Rod Style Motor Mount	12842-80
Flange Style Motor Mount	12842-90



1430 Northwest Boulevard
Vineland, NJ 08360

sales@aceglass.com
P: (800) 223-4524
F: (800) 543-6752