



# ACE-THRED™ FILTRATION APPARATUS

Using our famous Ace-Thred™ design, we are able to thread the filter funnel and filter flask together without the silicone stopper and troublesome clamping device. This design also eliminates leakage problems. Various porosities of filter discs are available. Fritted discs are easily removable, as it is not permanently sealed.

## 3700 Series Apparatus

- 500mL filter flask
- #7 Ace-Safe™ connector
- #25 to #25 Ace-Thred™ nylon adapter
- 25mm fritted disc
- 100mL graduated funnel

## 3702 Series Apparatus

- 1000mL filter flask
- #7 Ace-Safe™ connector
- #25 to #50 Ace-Thred™ nylon adapter
- 47mm fritted disc
- 500mL graduated funnel

## 3704 Series Apparatus

- Same great design as 3702 series
- Replaces the funnel with Instatherm® version
- Heats viscous materials for easier flow



### Also Available...

Contact us for all your filtration needs - both, large and bench scale applications.



## FILTRATION APPARATUS continued...

### 3708 Series Apparatus

- 2L filter flask
- #15 Ace-Safe™ connector
- #50 to #80 Ace-Thred™ PTFE adapter
- 75mm fritted disc
- 1L graduated funnel

### 3709 Series Apparatus

- 5L Duran bottle
- #15 Ace-Safe™ connector
- #80 Ace-Thred™ to GLS80 PTFE adapter
- 75mm fritted disc
- 1L graduated funnel

### 3710 Series Jacketed Funnels

- For 75mm membranes or filter paper
- Ideal for solvent recycling and other filter applications where temperature control is required
- 1L & 2L sizes available

**Contact us for...**

*Custom glass and repair services.*



1430 Northwest Boulevard  
Vineland, NJ 08360

sales@aceglass.com  
P: (800) 223-4524  
F: (800) 543-6752