



SERIES

6386 FILTER REACTORS

Jacketed Base Systems include:

- Jacketed Reactor Body w/O-Ring
(Capacities ranging from 150mL through 6000mL)
- Head and Clamp
- Bearing, Stir Shaft and Agitator
- Ace-Thred™ Bottom Adapter w/O-Ring
- Filter Support and Screen
- Bottom Drain Valve

Pressure Version adds:

- Condenser
- Addition Funnel
- Standard Taper joints replaced with Ace-Threds



JACKETED



Shown with optional Ace Impresario I, reactor control system.

Also Available...

Contact us for custom or large scale applications.



Contact your local sales rep today for more information

1-800-223-4524 | www.aceglass.com

SERIES

6384 FILTER REACTORS

Unjacketed Base Systems include:

- Unjacketed Reactor Body w/O-Ring
(Capacities ranging from 150mL through 6000mL*)
- Head and Clamp
- Bearing, Stir Shaft and Agitator
- Ace-Thred™ Bottom Adapter w/O-Ring
- Filter Support and Screen
- Bottom Drain Valve

Pressure Version adds:


- Condenser
- Addition Funnel
- Standard Taper joints replaced with Ace-Threds

Instatherm™ Filter Reactor Body

Same rugged, heavy wall reactor body, except Instatherm™ coated for safer, rapid heating. Operated with an Ace Temperature Controller. Available in the same capacities as the 6384 series unjacketed flask.



UNJACKETED



For safety, it is recommended that a pressure relief manifold be used on jacketed reactors.

Also Available...

Contact us for custom or large scale applications.



1430 Northwest Boulevard
Vineland, NJ 08360

sales@aceglass.com
P: (800) 223-4524
F: (800) 543-6752