



CELL CULTURE

FERNBACH FLASKS

Fernbach flasks allow maximal oxygen transfer to the culture medium, and provide a suitable surface via a large surface-to-volume ratio that enables increased cell divisions. Fernbach flasks also provide a hydrophilic environment, enabling optimal cell attachment.

3874 Series Flasks

- Beaded neck and accepts cotton plugs or rubber stoppers
- Triple baffled on the bottom edges to achieve maximal oxygen transfer to culture medium

3875 Series Flasks

- Long-style plain neck for stainless steel or polypropylene closures
- Triple baffled on the bottom to achieve maximal oxygen transfer to culture medium

3877 Series Flasks

- GL45 neck and accepts supplied polypropylene cap
- Triple baffled on the bottom edges to achieve maximal oxygen transfer to culture medium

3879 Series Flasks

- Long-style plain neck for stainless steel or polypropylene closures
- Flat bottom, no baffles

7622 Series Closures

- Various materials and colors available



Caution: Do NOT place on direct heat source or autoclave while capped.



Contact your local sales rep today for more information

1-800-223-4524 | www.aceglass.com