



JACKETED

WOULFF 3-NECK BOTTLE

The Woulff Jacketed Three Neck Flask is designed to provide scientists, researchers, and professionals in the chemical industry with a superior platform for conducting complex reactions. Its jacket ensures temperature control, enabling users to carry out reactions that demand regulation of heat.

Key Features

- Enhanced temperature control with the use of the integrated jacket.
- Versatile configuration through three necks.
- Durable build as you have come to expect from Ace Glass's high-quality materials.
- Safety first by minimizing overheating and unwanted reactions due to the jacketed design.
- Integrated Ace-Thred™ jacket #11 ports supplied with Ace-Safe connectors for 3/8" ID tubing.

Various Uses

- Intermediary flask for vacuum work to help collect solvents and/or control vacuum.
- They may also be used for washing gases or dissolving them in solvents.
- In a rotovap system, a Woulff bottle can be used to prevent solvents from damaging the pump.

Also Available...

Non-Jacketed Woulff Heavy Wall Bottle

