



# STREAMFLOW™ SOLUTION TRANSFER BOTTLE

Upgrade your laboratory operations with our StreamFlow™ Solution Transfer Bottle, available in a range of sizes. This advanced storage solution features a hose connection and a vented cap with an expanded PTFE membrane, enabling effortless and efficient transfer of media while allowing for controlled gas exchange. It is equipped with a drip-free pouring ring and the expanded PTFE membrane in the vented cap facilitates sterile gas exchange while effectively blocking high surface tension liquids.

## Effortless Solution Transfer:

- Seamless media transfer with the hose connection
- Expanded PTFE membrane in the vented cap enables controlled gas exchange
- Avoid spills and contamination with the drip-free pouring ring

## Controlled Gas Exchange:

- Expanded PTFE membrane facilitates sterile gas exchange
- Ensures protection against contaminants while allowing for necessary gas exchange
- Ideal for sensitive cultures or samples requiring controlled atmospheric conditions

## Reliable Performance:

- Autoclavable properties for easy sterilization
- Durable construction ensures long-lasting use
- Suitable for various laboratory settings, including research, pharmaceutical, and biotech industries

## Range of Sizes:

- Choose the perfect capacity for your specific needs

