



# ***FREEZE DRYING*** **FLASK w/SNAP-ON CAP**

Complete “fast-freeze” flask with snap-on rubber cap and plain or stopcock PTFE adapter. Flasks are designed for ease of handling and faster loading. Wide mouth provides easy access to inner walls. Snap-on black rubber cap is grease free; one size fits all capacity flasks. Plain straight adapter for connecting flask cap to lyophilizer manifold is removable.

## ***7035 Series Flasks***

- Compatible with all major brands of freeze dryers and manifolds
- Uniform wall weight for good temperature transfer
- Wide mouth provides ease in getting product in and out of the flask
- Black rubber cap matches O.D. and I.D. on the flask
- PTFE adapter tubes available plain or with stopcock
- Caps fit standard lyophilizer manifolds



***Contact us for...***

*Custom glass and repair services.*



Contact your local sales rep today for more information

1-800-223-4524 | [www.aceglass.com](http://www.aceglass.com)

# ***FREEZE DRYING*** **FLASK w/TWISTLOR SEAL**

These are just a few of the more widely used glassware components to replace your existing glass. All our rotary evaporator replacement glass is manufactured in the USA, to exacting specifications and standards.

## ***7030 Series Flasks***

- O-Ring for vacuum-tight seal
- Dual purpose: Rotary evaporator and freeze drying flask
- Adapter heads will fit all capacity flasks
- Adapter head comes with nylon retainers and FETFE O-Ring
- Note: These flasks are NOT for use in elevated temperatures



### ***Also Available...***

Contact us for all your freeze drying and rotavap glassware needs, for both large and bench scale applications.



1430 Northwest Boulevard  
Vineland, NJ 08360

[sales@aceglass.com](mailto:sales@aceglass.com)  
P: (800) 223-4524  
F: (800) 543-6752