



FREEZE DRYING

FLASK w/SNAP-ON CAP

Complete “fast-freeze” flask with snap-on rubber cap and plain or stopcock PTFE adapter. Flasks are designed for ease of handling and faster loading. Wide mouth provides easy access to inner walls. Snap-on black rubber cap is grease free; one size fits all capacity flasks. Plain straight adapter for connecting flask cap to lyophilizer manifold is removable.

7035 Series Flasks

- Compatible with all major brands of freeze dryers and manifolds
- Uniform wall weight for good temperature transfer
- Wide mouth provides ease in getting product in and out of the flask
- Black rubber cap matches O.D. and I.D. on the flask
- PTFE adapter tubes available plain or with stopcock
- Caps fit standard lyophilizer manifolds



Contact us for...

Custom glass and repair services.



Contact your local sales rep today for more information

1-800-223-4524 | www.aceglass.com

FREEZE DRYING **FLASK w/TWISTLOR SEAL**

These are just a few of the more widely used glassware components to replace your existing glass. All our rotary evaporator replacement glass is manufactured in the USA, to exacting specifications and standards.

7030 Series Flasks

- O-Ring for vacuum-tight seal
- Dual purpose: Rotary evaporator and freeze drying flask
- Adapter heads will fit all capacity flasks
- Adapter head comes with nylon retainers and FETFE O-Ring
- Note: These flasks are NOT for use in elevated temperatures



Also Available...

Contact us for all your freeze drying and rotavap glassware needs, for both large and bench scale applications.



1430 Northwest Boulevard
Vineland, NJ 08360

sales@aceglass.com
P: (800) 223-4524
F: (800) 543-6752